

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 1 of 15 |

1. Purpose:

Assure reliability of handheld and cable on the Gen 1B Blood Flow device going to clinical Trial. ASM, Tower Gen 1B 7000-0094

2. Scope:

This document applies to part 7000-0089 ASM, wand 1B bundle and 7000-0122 ASM, Wand 1B

3. Background (optional):

According to IEC 60601 handheld devices and cable must be robust to expected usage.

4. Reference Documents:

Drop test from IEC 60601-1 drop test for handheld equipment.

5. Testing

1. Drop test

15.3.4.1 HAND-HELD ME EQUIPMENT

HAND-HELD ME EQUIPMENT, ACCESSORIES and ME EQUIPMENT parts that are HAND-HELD shall not result in an unacceptable RISK as a result of a free fall.

Compliance is checked by the following test.

The sample to be tested, with any SAFE WORKING LOAD in place, is allowed to fall freely once from each of three different starting orientations encountered during NORMAL USE from the height at which the ME EQUIPMENT, ACCESSORY or ME EQUIPMENT part is used (as specified in the ACCOMPANYING DOCUMENTS), or

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 2 of 15 |

from a height of 1 m, whichever is greater, onto a 50 mm ± 5 mm thick hardwood board (hardwood > 600 kg/m³) lying flat on a concrete or a similar rigid base.

After the test, the HAND-HELD ME EQUIPMENT and, ACCESSORY or ME EQUIPMENT parts shall not result in an unacceptable RISK.



Drop 1 distance

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 3 of 15 |



dropped 3 different orientations 3 times each.
9 total drops from 1 meter

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 4 of 15 |



no physical damage

2. rollover crush test

weigh one wheel of a tower or cart, both are the same weight.

rollover cable 11 times with a representative tower or cart.

verify function of cable before and after test by making a blood flow measurement

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 5 of 15 |



| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 6 of 15 |



No physical damage after the crush test and drop test.



laser beam ok after crush and drop test

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|---------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 7 of 15 |



all 4 cameras working properly

measurement working properly- OK

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|



ENGINEERING STUDY REPORT

| | |
|-----------------------------------|---------------|
| Doc #: | D00036 |
| Rev: | 01 |
| Reliability test, Gen 1B handheld | Page: 8 of 15 |

3. Pull test

From IEC 60601

The cord is subjected 25 times to a pull on the sheath of the value shown in Table 18. The pulls are applied in the most unfavourable direction without jerks, each time for 1 s.

Immediately afterwards, the cord is subjected for 1 min to a torque of the value shown in Table 18.

Table 18 – Testing of cord anchorages

| Mass (<i>m</i>) of ME EQUIPMENT kg | Pull N | Torque Nm |
|---|-----------|--------------|
| $m \leq 1$ | 30 | 0,1 |
| $1 < m \leq 4$ | 60 | 0,25 |
| $m > 4$ | 100 | 0,35 |

A cord anchorage that allows the cord sheath to be longitudinally displaced by more than 2 mm or the conductor ends to move over a distance of more than 1 mm from their normally connected position is considered to fail.


CREEPAGE DISTANCES and AIR CLEARANCES that are reduced below the values specified in 8.9 constitutes a failure.

Attempt to push the cord into the ME EQUIPMENT or the MAINS CONNECTOR. If the cord can be pushed into the ME EQUIPMENT or the MAINS CONNECTOR to such an extent that the cord or internal parts are damaged, the cord anchorage is considered to fail.

**ENGINEERING STUDY REPORT**

| | |
|-----------------------------------|---------------|
| Doc #: | D00036 |
| Rev: | 01 |
| Reliability test, Gen 1B handheld | Page: 9 of 15 |

Anchor picture bulkhead side
cross section of model shown below:



| | | | |
|---|-----------------------------------|--------|----------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 10 of 15 |




102.45N on tower bulkhead,
Passed

Anchor picture handheld side

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|----------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 11 of 15 |



39.6 N on wand handle, pull test ok.

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|----------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 12 of 15 |

4. cleaning wipe test-
 pull wipe from CaviWipe container and wipe the entire wand surface and length 4 times.
 CaviWipe recommended by HH IRB

6. Study Summary:

Summary of tests made to the wand and cable

1. drop test hand held

An hand held 7000-0122 dropped from 1 meter onto a wooden board $\frac{3}{8}$ thick on the concrete floor 9 times total (3 times, 3 orientations).- no cracking

2. cable roll over

The cable was run over by the cart 11 times- no degradation, visual or functional

3. Pull test

a. pull test in bulk head to tower 3x
 passed, 100N ok

b. pull test on handheld 3 x
 passed, 39.6 N ok

4. wipe down/ cleaning

CaviWipe, 4 wipes handle and cable.
 no degradation seen, wipedown ok

The silicone sheathing and ABS handle were wiped with CaviWipe 4 times - no degradation visual or functional.

6.1. Equipment/Materials

6.1.1. Taylor scale lbs.

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|----------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 13 of 15 |

6.1.2. Chatillon force gauge DFS2-100 WS2 35955

6.1.3. Gen 1A cart

7. Results and Analysis:

| test # | Test | Result |
|--------|------------|--------|
| 1 | Drop test | Pass |
| 2 | Crush test | Pass |
| 3 | Pull test | Pass |
| 4 | Wipe test | Pass |

8. Conclusion/Recommendation:

Summary: The wand and cable is low risk for failing foreseen usage

The wand 1B is a passive device, except for the LED and tactile switch.

The tactile dome switch is a very robust design used in millions of consumer and medical devices.

The LED is a very robust solid-state device. The fibers are passive devices; any failure will cause less light to get through the laser.

9. Appendix List: *Chatillon calibration certification*

Appendix 1

| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|

| | | | |
|---|-----------------------------------|--------|----------|
|  | ENGINEERING STUDY REPORT | Doc #: | D00036 |
| | | Rev: | 01 |
| | Reliability test, Gen 1B handheld | Page: | 14 of 15 |




| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|



ENGINEERING STUDY REPORT

| | |
|-----------------------------------|----------------|
| Doc #: | D00036 |
| Rev: | 01 |
| Reliability test, Gen 1B handheld | Page: 15 of 15 |



| | |
|------------------------|--|
| TMPL #: D0002 Rev A | Printed copies are considered UNCONTROLLED documents after the date printed. Printed on October 8, 2020.CONFIDENTIAL AND PROPRIETARY - Openwater INC. |
|------------------------|--|